

EMV+[®]

CLEANING INSTRUCTIONS



Process for Cleaning the EMV+

- Keep the ventilator and its accessories clean at all times. Never allow grease or oil to enter the system or coat its components.
- Clean the device at regular intervals and maintain up-to-date records of inspections, cleaning, and maintenance.
- Take care to prevent liquids from entering the ventilator. Never submerge the ventilator and avoid using excessive amounts of water that might enter the device. Dry all exposed parts following use in wet environments.
- Clean the device's housing and hose connections with a damp, soapy cloth.
- For general decontamination and cleaning, apply a 10% bleach solution with a damp cloth.
- Do not clean the device with abrasives or chlorinated hydrocarbon cleansers, which damage the housing and interface lens.
- After cleaning, thoroughly dry the device with a lint-free cloth. Make sure that all exposed surfaces are cleaned and dried.
- Refer to the Operator's Guide for full cleaning and operating instructions.

WARNING!

Never use oil or grease of any kind with oxygen or compressed gas equipment.

Post-Contaminated Environment Cleaning

If you have used the ventilator in an environment where it may have been exposed to contamination from a hazardous materials accident, mass epidemic, or weapon of mass destruction, we recommend that you follow these guidelines:

- Always follow the decontamination procedures specified by the local Incident Command Safety Officer.
- You should clean and decontaminate the equipment as soon as possible after use. Personnel should always wear the appropriate Personal Protective Equipment while decontaminating equipment.
- Review to the aforementioned cleaning instructions.
- Cleaning and sterilizing practices may vary between institutions. We suggest that each facility have in place a procedure for the cleaning and disinfection of its medical equipment and that these procedures to be consulted for further guidance.