

# Reduce the Risk of Cross-Contamination with Disposable Accessories

Note: This list is for international use only. Some products are not available in the United States. U.S. Military/Government customers, please contact your local ZOLL® representative for a list of disposable accessories available in the U.S.

To reduce the risk of cross-contamination, ZOLL offers a variety of disposable products for use with:

## Monitors/Defibrillators

- Propaq® M and MD
- X Series®



## Ventilators

- EMV+®
- Z Vent®



The use of disposable accessories can help to reduce:

- risk of indirect cross-contamination
- risk of direct contamination during the cleaning process
- often lengthy reprocessing procedures
- use of valuable disinfectants

## DISPOSABLE MONITOR/DEFIBRILLATOR ACCESSORIES

### BLOOD PRESSURE CUFFS

- Every cuff size, from small infant through thigh cuffs
- Color-coded for fast, accurate size selection
- Rotatable port improves patient comfort and cuff durability
- Meets all the latest clinical guidelines for proper fit from AAMI and AHA
- Durable enough for extended patient stays — rated to over 500 life cycles (includes on/off arm placements and inflations)
- Latex-free



SOFT-07-2MQ	Disposable Cuff, soft infant double tube (box of 20)
SOFT-08-2MQ	Disposable Cuff, soft sm child double tube (box of 20)
SOFT-09-2MQ	Disposable Cuff, soft child double tube (box of 20)
SOFT-10-2MQ	Disposable Cuff, soft sm adult double tube (box of 20)
SOFT-11-2MQ	Disposable Cuff, soft adult double tube (box of 20)
SOFT-11L-2MQ	Disposable Cuff, soft adult long double tube (box of 20)
SOFT-12-2MQ	Disposable Cuff, soft lg adult double tube (box of 20)
SOFT-12L-2MQ	Disposable Cuff, soft lg adult long double tube (box of 20)
SOFT-13-2MQ	Disposable Cuff, soft thigh double tube (box of 20)

## SPO<sub>2</sub> SENSORS


### Two-piece M-LNCS – Rainbow® Technology

- Single-patient use sensors have a flat cable with smooth edges that lie comfortably on a patient's hand
- Appropriate for short- or long-term monitoring of adult, pediatric, infant, and neonatal patients
- Designed for finger, toe, thumb, hand, or foot application, offering secure positioning of optical components
- Minimize the risk of cross-contamination



## SPO<sub>2</sub> ONLY

(For use with 8000-0331 or 8000-0330 reusable patient cable)



**M-LNCS Adhesive Sensor**  
8000-000475  
(ADTX Adult, Box of 20)  
8000-000476  
(PDTX-3 Pediatric, Box of 20)



**M-LNCS Inf-3 Adhesive Sensor**  
8000-000477 (Infant, Box of 20)



**M-LNCS Neonatal Preterm Adhesive Sensor**  
8000-000478  
(NeoPt-3 Neo, Box of 20)

## SPO<sub>2</sub>, SPCO, SPMET

(For use with 8000-0331 or 8000-0330 reusable patient cable)



**Rainbow Single-Use Sensors**  
8000-0336 (Patients > 30kg)  
8000-0339 (Patients 10-50kg)



**Rainbow Single-Use Sensors**  
8000-0337  
(Patients < 3kg, > 30 kg)  
8000-0340  
(Infant Patients 3-30kg)



**R1-25 Butterfly Adult Adhesive Sensors**  
8000-000468  
SpHb, SpMet, SpO<sub>2</sub>  
(Box of 10)



**R1-25L Adult Adhesive Sensors**  
8000-000462  
SpHb, SpMet, SpO<sub>2</sub>  
(Box of 10)



**R1-20 Butterfly Pediatric Adhesive Sensors**  
8000-000469  
SpHb, SpMet, SpO<sub>2</sub>  
(Box of 10)



**R1-20L Infant Adhesive Sensors**  
8000-001463  
SpHb, SpMet, SpO<sub>2</sub>  
(Box of 10)

## TEMPERATURE SENSORS

ZOLL offers single-use temperature probes for use in different clinical applications



**YSI Single-Use Adult Esophageal Rectal Probe**  
8000-0672



**YSI Single-Use Adult Skin Temperature Probe**  
8000-0673



**Disposable Temp Sensor Adapter Cable**  
8000-0674

## EtCO<sub>2</sub> SAMPLING LINES

ZOLL offers a full array of EtCO<sub>2</sub> sampling lines that meet every clinical need for intubated and non-intubated patients.

FilterLine non-intubated sampling lines feature oral and nasal sampling as well as an innovative supplemental oxygen delivery system. All FilterLine non-intubated sampling lines can be used with bi-level ventilation and CPAP.

- Integrated, single-piece airway gas sampling lines include the patient cannula or airway adapter, water trap, and a sterilizing grade filter that is designed to reduce the risk of biohazard contamination of the patient monitor
- Uni-junction technology is designed to ensure proper sampling for patients who alternate between oral and nasal breathing
- Oral scoop design has wide surface area to provide breath capture even in the presence of shallow breathing
- Oxygen flows from holes in front of nasal prongs and improves humidity of O<sub>2</sub>
- Delivers O<sub>2</sub> flow up to 5 L/min and can deliver above 5 L/min using external O<sub>2</sub> delivery, e.g. through a mask, without affecting CO<sub>2</sub> measurement values
- Broad portfolio of sampling lines available for various clinical situations
- One size fits adolescents through adults; for some types of sampling lines, dedicated pediatric products are also available



**Oridion Filderline Set Adult/Pediatric**  
8300-0520-01 (Case of 25)



**Oridion Vitaline H Set Adult/Pediatric**  
8300-0523-01 (Case of 25)



**Oridion Filterline H**  
8300-0521-01 (Ad/Ped, Set of 25)  
8300-0522-01 (Inf/Neo, Set of 25)



**Oridion Smart CAPNO line Plus O<sub>2</sub>**  
8300-0524-01 (Adult, Case of 25)  
8300-0525-01 (Pediatric, Case of 25)

## DISPOSABLE VENTILATOR ACCESSORIES

### HEAT MOISTURE EXCHANGERS

Heat moist exchanger with filter; provides humidification to the patient during ventilation.  
Not for use with CPAP or BL.



#### HME Adult

465-0025-00 (Case of 20)

820-0053-00 (Case of 50)



#### HME Pediatric Dead Space <25ML

820-0109-00 (Case of 20)



#### HME Adult Dead Space <75ML

820-0108-00 (Case of 20)



#### HME Infant Dead Space <10ML

820-0110-00 (Case of 30)



#### Adult Circuits

6-foot disposable, single use

820-0106-15 (Case of 15)



#### Bacterial/Viral Filter

May be used on the fresh gas intake

465-0024-00 (Case of 50)

### PULSE OXIMETER SENSORS (LNCS)

- Appropriate for short- or long-term monitoring of adult, pediatric, infant, and neonatal patients
- Designed for finger, toe, thumb, hand, or foot application, offering secure positioning of the optical components
- SoftTouch sensors are designed for patients with fragile skin and feature a non-adhesive wrap made from soft foam material, with hook-and-loop attachment



#### LNCS Reusable SpO<sub>2</sub> Patient Cable

8000-0298 (4 feet)

8000-0293 (10 feet)



#### LNCS Neonatal SpO<sub>2</sub> Adhesive Sensor

8000-0324 (Neonatal, Box of 20)



#### LNCS Disposable SpO<sub>2</sub> Sensors

8000-0320 (Adult, Box of 20)

8000-0321 (Pediatric, Box of 20)



#### LNCS SpO<sub>2</sub> Adhesive Sensors

8000-0322 (Infant, Box of 20)

8000-0323 (Neonatal, Box of 20)