



Patient-Simulation System



Our patient-simulation systems enables you to run **more realistic simulation trainings**. Because of the consequent design, real-life functionalities and ease of use our systems support you to achieve a **higher quality training results**. You and your students will have fun and you will love it.

Our **patienten-simulation systems** consist of

- an iPad Pro ECG-monitor,
- an instructor-controller for iPad,
- an assessment-application for iPod and iPhone as well as
- a case including a set of cable.

qube3

monitor interface of **corpuls3**

qube15

monitor interface of **LIFEPAK15**

qube7

monitor interface of **DEFIGARD® Touch7**

qubeX

world famous monitor interface

dynamic **12-Lead ECG with possibility to print**

defibrillation (AED/manual) with feedback

pacemaker / pacemaker

CPR-mode with CPR-fraction measurement, metronome, recording no-flow-time, etc.

instructor control for iPad

assessment-application for iPod and iPhone

in cooperation with **GS Elektromedizinische Geräte G. Stempel GmbH**

in cooperation with **Stryker GmbH & Co. KG**

in cooperation with **SCHILLER**

—



Full functional **ECG-monitor interface**, e.g. with CPR-mode, defibrillation and pacemaker



Auscultation functionality, blood glucometer, thermometer, checklists, SOP, etc. via the assessment application



Seamless integration with our **qubeAV** mobile debriefing system - turn this flyer around



Browser-based scenarioeditor, image and video library, control of vital parameters by instructor



Checklists and SOP available for **documentation** of simulation training

Debriefing-System



Debriefing is one of the most relevant part of simulation training. With our **qubeAV** debriefing systems you will be able to support the students during and run the debriefing after the training.



Mobile debriefing-system with audio-video-recording for simulation training



3 camera channels + 1 channel to intergate qube³ and qube¹⁵ patient simulation system, Student view - live and on demand



Recording the training allows analysis of the training scenario from all angles



Multiple marker-sets defined by your own, controlled via the marker keyboard, track pad or iPad



Due to the **mobility** of the system you are totally **flexible** and you can **choose** if you like to train in the classroom, on the ward, in the ambulance or outside.



Seamless integration with our **qube³ & qube¹⁵** patient simulation system - turn this flyer around

qubeAVpro

qubeAVgo

professional debriefing-system

mobile debriefing system

macbook based software

iPad based Software Application

3 HD-cameras

360° camera

stream external monitors

stream SKILLQUBE monitors

Voice of God

smart Microphone

Marker Keyboard

set Marker with ApplePen



Medical Simulator
INNOVACIÓN EN EDUCACIÓN

Medical Simulator Spain
Ctra de Pozuelo a Majadahonda Km. 1800
28223 Pozuelo de Alarcón | Madrid | Spain
T 902 196 788
T (+34) 91 382 08 88
F (+34) 91 381 98 80

www.medical-simulator.com
info@medical-simulator.com
@Simulacion_Medi
medical-simulator
medicalsimulator
@medical_simulator